

BetterLiving

Non-Slip Bedside Safety Crash Mat



helps to prevent patient injury due to falls from bed

sectioned and hinged for easy fold and storage

includes slip resistant textured underside to minimise slip hazard

reflective strip for easy night-time visibility

This Safety Crash Mat can be placed on either side of the bed as a protective cushion to help prevent injury for patients who roll out of bed. Especially useful where side rails cannot be used.

The new non-slip feature helps to minimise the risk of injury should the mat be inadvertently stood on by patient or carer. The unique textured underside prevents the mat slipping when downward pressure is applied however still allows the mat to be repositioned with the foot for easy manoeuvrability.

The Crash Mat includes a reflective strip to minimise trip hazards in low light visibility. Comprises three sections of safety foam covered by hinged nylon for easy folding and storage.

Code	BSBA7272
Thickness	50 mm
Width	600 mm
Length	1750 mm
Artg	155475
Warranty	2 year



HospEquip Australia | www.hospequip.com.au | EMAIL info@hospequip.com.au

WESTERN AUSTRALIA
21 Tulloch Way, Canning Vale WA 6155
T: (08) 9456 1661

VICTORIA
2 Bate Dr, Braeside VIC 3195
T: (03) 9580 1055

SOUTH AUSTRALIA
62 Leader St, Forestville SA 5035
T: (08) 7225 8720

QUEENSLAND
5/14 Ashtan Pl, Banyo QLD 4014
T: (07) 3181 5971