

LINAK JUMBO BATTERY PACK

**SPECIAL
\$180 +
GST**

The BAJ is a lead acid battery for use together with the modular JUMBO™ patient lift system. It offers variable and reliable battery capacity together with easy mounting, charging and replacement.

The BAJ Linak battery pack will work with most lifting hoists and standing machines from Direct Healthcare Group, System RoMedic, Bure, Allegro, Kerry, Invacare, Pro-Med, Liko, Haycomp, Novis, Liftaid, Oxford & Aspire.

We recommend replacing the battery every 5 years for optimum battery life and performance.

CHARGING

BAJ can be charged by:

- Charger CHJ2
- Control box CBJ2, CBJ Care, COBO
- BAJ with integrated DC plug can also be charged by use of the eternal charger (CH01)



**ORDER CODE
BAJ1000001**

Duty Cycle	Ambient Temperature	Storage Temperature	Relative Humidity	Atmospheric Pressure	Height Above Sea Level	Fammability Rating	Approvals
10% or 2 min. continuous use then 18 min. not in use	+5 to +40° C	-15° to +40° C	20% to 80% - non-condensing	700 to 1060 hPa	Max. 3000 meters	UL94-V0	IEC60601-1 ANSI/AAMI ES60601 CAN/CSA-22.2 No 60601-1

SPECIAL PRICE EXCLUDES DELIVERY & TAXES

CONTACT US

WA 📞 (08) 9456 1661
✉️ saleswa@hospequip.com.au

VIC 📞 (03) 9580 1055
✉️ vicadmin@hospequip.com.au

🌐 hospequip.com.au

SA 📞 (08) 7225 8720
✉️ sareception@hospequip.com.au

QLD 📞 (07) 3448 9243
✉️ qldadmin@hospequip.com.au

NSW 📞 1300 467 737
✉️ nswadmin@hospequip.com.au